

Novel Method for High Molecular Weight (HMW) DNA Extraction


P. Giron Koetsier, Barbara Taron, Luo Sun, Danielle Freedman, Eileen Dimalanta and Eric Cantor | New England Biolabs, Inc.



INTRODUCTION


Extraction of high molecular weight (HMW) genomic DNA (gDNA) is currently a pain point for many researchers, with existing solutions taking hours, or even days, to complete. Under the Monarch® Nucleic Acid Purification product portfolio, we have developed a rapid and reliable process for extracting intact HMW genomic DNA from mammalian tissues, cultured cells, blood, and bacteria. Utilizing an optimized process that combines gentle lysis with tunable fragment length generation, followed by precipitation of extracted DNA onto the surface of glass beads, the prep proceeds rapidly and utilizes standard laboratory equipment. DNA size ranges from 50-250 kb for the standard protocol and can be adjusted to produce longer DNA up to >1 Mb. Purified DNA is recovered in high yield with excellent purity, including nearly complete removal of RNA. Prep time is approximately 30-90 minutes, depending on sample type. The purified HMW gDNA is suitable for a variety of downstream applications including long-read sequencing and genome mapping.

WORKFLOW: CELLS & BLOOD



PERFORMANCE: TUNABLE LENGTH, YIELD & PURITY

Use of varying agitation speeds during lysis produces controlled fragment length of extracted HMW genomic DNA from cells, blood and soft organ tissue.



Preps were performed on aliquots of 10⁶ cells (HEK293, duplicates), 500 µl human blood (duplicates), or mouse kidney (10 mg, single). Samples were agitated at the indicated speed during the lysis step to control the fragmentation of the DNA. Equal amounts of DNA from the replicates (cells: 500 ng; kidney 300 ng) were resolved by PFGE. Yield and purity of the individual preps are shown in the accompanying tables.

ONT SEQUENCING METRICS*

Table 1: Cells & Blood

	HEK293		Human blood		K562 cells	Chicken blood	Pig blood
	Single run	Single run	Single run	Single run	Multiplex	Multiplex	Multiplex
Read length N50 (kb):	45.4	40.4	51.4	46.5	41	37.9	37.6
Total bases (Gb):	8.1	12.2	8.1	11.6	0.37	0.41	0.41
Mean read length (kb):	21.3	19.2	24.7	21.5	16.9	16	16.6
Median read length (kb):	10.4	9.7	12.6	10.1	7.1	7	7.9
Mean read quality:	12.8	13.2	13.3	13.4	13.2	12.9	13
Median read quality:	13.2	13.7	13.8	13.9	13.6	13.4	13.5

Table 2: Tissues & Bacteria

	Kidney	Kidney	Mouse brain	Mouse liver	Mouse muscle	B. cereus	E. coli
	Single run	Single run	Multiplex	Multiplex	Multiplex	Multiplex	Multiplex
Read length N50 (kb):	33	44.6	40.9	42.3	37.9	43.4	48
Total bases (Gb):	7.1	4.5	0.12	0.36	0.13	0.23	0.31
Mean read length (kb):	15.9	27.1	21.9	27.4	18.9	22.2	29.4
Median read length (kb):	8.8	23.2	14.8	23.7	10.9	13.6	23.8
Mean read quality:	13	13	12.9	13.2	12.8	12.5	12.7
Median read quality:	13.4	13.5	13.3	13.7	13.4	13.1	13.2

Table 3: Comparison of Monarch to Circulomics® Nanobind®

*Identical samples in multiplex runs	Fresh K562 Cells (7x10 ⁵)						Fresh HEK293 Cells (2x10 ⁶)						Human Blood (200 µl)					
	Monarch			Circulomics			Monarch			Circulomics			Monarch			Circulomics		
	1	2	3	1	2	3	1	2	3	1	2	3	1	2	3	1	2	3
Read length N50 (kb):	41	40.2	40.8	37.6	41.7	39.8	40.7	41.8	40	40.8	41.5	40.7	38.1	37.8	38.3	32.9	31.6	30.7
Total bases (Gb):	0.37	0.77	0.19	0.47	0.43	0.29	0.42	0.44	0.36	0.31	0.35	0.21	0.43	0.46	0.28	0.47	0.36	0.44
Mean read length (kb):	16.9	16.4	17.4	12.5	16	13.6	18.4	18.6	17.8	18.4	18.2	17.8	17.1	16.7	16.7	16	15.7	15.5
Median read length (kb):	7.1	6.9	7.4	4.5	6.2	4.8	8.6	8	8.1	8.9	8.5	8.2	7.9	7.6	7.5	8.3	8.3	8.2
Mean read quality:	13.2	13.1	13.1	13.2	13.2	13.2	13.3	13.3	13.3	13.1	13	13	13.3	13.3	13.3	13.3	13.3	13.3
Median read quality:	13.6	13.5	13.5	13.5	13.5	13.5	13.7	13.7	13.8	13.6	13.5	13.5	13.7	13.7	13.8	13.7	13.7	13.7
Yield (µg)	8.1	8.2	8.2	4.9	4.9	5.2	20.5	22.9	23.4	15.9	15.0	17.1	11.4	10.8	11.6	11.6	10.4	12.4
260/280	1.85	1.87	1.87	1.86	1.83	1.72	1.86	1.85	1.86	1.84	1.84	1.84	1.85	1.87	1.87	1.81	1.80	1.85
260/230	2.43	2.52	2.31	2.24	2.55	1.98	2.43	2.54	2.48	2.28	2.01	2.34	2.43	2.49	2.49	1.85	1.95	1.86

*LSK109, FLO-MIN106D; load 800-1000 ng HMW DNA per flow cell for optimal molar ratios

HIGHLIGHTS

- Rapid protocols with minimal hands-on time
- Maximal yield and purity
- Purified DNA ranges from 50 kb to > 1 Mb
- DNA size from cells, blood, bacteria and soft organ tissues can be tuned based on agitation speeds
- Standard protocols yield DNA 50-250 kb for Oxford Nanopore Technologies® (ONT) sequencing
- Highly reproducible
- Low RNA content
- RNase A, Proteinase K, microtube pestles included


Fast workflows:

- Cells: ~30 min
- Blood: ~60 min
- Tissue: ~90 min
- Bacteria: ~90 min

High N50s:

- Cells: 40-55 kb
- Blood: 40-55 kb
- Tissue: 35-45 kb
- Bacteria: 40-48 kb

WORKFLOW: TISSUE & BACTERIA



PERFORMANCE: COMPARISON HMW DNA EXTRACTION KITS

Monarch HMW DNA Extraction Kit produces high molecular weight DNA with excellent yields and purity as compared to other commercially available kits

