



Validation and quality control of a molecular cloning experiment using *de novo* assembly of Oxford Nanopore reads

wf-clone-validation — for rapid verification of artificial constructs such as plasmids, bacterial artificial constructs (BACs), and fosmids

Contact: support@nanoporetech.com. More information can be found at: labs.epi2me.io/workflows/wf-clone-validation/. Data used in this analysis is available to download from: epi2me.nanoporetech.com/plasmid_2025.04/




Fig. 1 a) Legacy methods vs clone validation end-to-end workflow b) sample and library preparation diagram c) sequencing and basecalling devices d) EPI2ME™ application for analysis.

End-to-end solution

Legacy methods for validation of cloned vectors usually only consider a small section of the construct. Therefore, they can miss changes in the vector backbone, such as recombination and mutations, or insertion failure (Fig. 1a). The end-to-end clone-validation workflow overcomes these issues by assembling whole constructs, ranging in size from 1,500 to 300,000 bp (Fig. 1b, c, d). To demonstrate its utility, we have used a 96 sample plasmid dataset, which contains two sample prep replicates and two technical replicates, across four MinION Flow Cells. The samples were prepared following the the Rapid sequencing V14 — Plasmid sequencing protocol (SQK-RBK114.96) and sequenced on a MinION™. Basecalling was performed using the Dorado dna_r10.4.1_e8.2_400bps_hac@v5.0.0 HAC model and the FASTQ files were analysed with the EPI2ME wf-clone-validation workflow (v1.8.0) on a standard laptop (16 GB RAM; 8 CPUs).




Fig. 2 a) Running analysis in the EPI2ME application b) sample status table indicating if assembly was successful, assembly length, mean quality, and if the insert and construct recovered were expected compared to the reference.

Comprehensive assemblies

The sequencing output files were directly input into the wf-clone-validation workflow in the EPI2ME application. A sample sheet with per-sample approximate construct sizes, references, and insert references was provided for additional quality control after the *de novo* assembly. The workflow completed in less than an hour and progress was monitored in the workflow overview tab. The workflow produced a per-sample construct FASTA, insert FASTA, plasmid annotations in BED and GBK formats, and a report. The sample status table in the report (Fig. 2b) showed that all assemblies completed successfully, and we can see the read lengths and mean qualities. However, for sample_11, although a plasmid was assembled, when compared to the provided reference it did not meet the acceptance criteria (default 95% coverage, 99% identity) so it showed as failed.




Fig. 3 a) pLannotate annotation plot b) pLannotate feature table c) raw and downsampled read statistics.

Detailed annotations

As can be seen for the example of sample_11, the pLannotate annotation plot in the report provides a visual representation of each construct, confirming whether the expected genes or features are present and in the correct orientation relative to the plasmid (Fig. 3a). If a reference was provided, the orientation of the plasmid will match that of the reference. A full feature table (Fig. 3b) offers additional detail, including percent identity, match length, and descriptive annotations. The read summary section (Fig. 3c) provides insight into read quality and length distributions; for instance, in sample_11, a sharp peak at approximately 7 kb indicates that most reads span the full plasmid. The downsampled data used for assembly can be adjusted by modifying the coverage and minimum quality score settings. Sample_11 has good read lengths and qualities indicating no issues with the sequencing experiment.




Fig. 4 a) Assembly QC from b) multiple sequence alignment c) self-alignment dot plot.

Quality control

The QC table in the report (Fig. 4a), which compares the assembled constructs with the provided reference, shows that sample_11 only has 30% of coverage, which is less than the expected ≥95%. A multiple sequence alignment between provided insert references and the sample confirms the sample_11 insert is not as expected (Fig. 4b). The matching coverage is likely because the vector matches but the insert is not as expected. There are various other QC metrics available to help inspect constructs, including a self-alignment dot plot: repetitive regions appear as diagonal lines off the main diagonal, indicating internal duplications or repeats. These can be used to check repetitive regions have assembled as expected (Fig. 4c). To learn more and try out the workflow, visit: epi2me.nanoporetech.com/epi2me-docs/workflows/wf-clone-validation/.