Min**ION** Mk 1B

Quick start guide

This Quick start guide contains everything you need to set up your MinION™ Mk1B and to check that the device is ready for use.

For detailed information and troubleshooting, view the user manual.

Pre-installation

Before using the device, familiarise yourself with the following:

IT requirements

community.nanoporetech.com/to/minion-it

MinION Mk1B user manual

community.nanopore.com/to/minion

Safety and regulatory information

community.nanoporetech.com/to/safety

What's in the box:

MinION Mk1B

CTC

Configuration Test Cell

USB Type-A cable*

*Requires USB 3.0 speeds or greater

1

Overview

To start, check that your hardware is working.

Install
MinKNOW™

Connect the
MinION to your
computer

Perform a
hardware
check


Connect with
the Nanopore
Community

2

Install and log in to MinKNOW

Computer:

1. Check that the computer you are using meets the IT requirements.
2. Power on your computer.
3. Install the latest **MinKNOW software for the MinION Mk1B and the PromethION 2 Solo**

devices from the Software Downloads page on your computer. Once the software is installed, click the Nanopore wheel icon that appears on the desktop.

4. Log in to MinKNOW using your Nanopore account details.


Software downloads
community.
nanoporetech.
com/downloads

3

Set up your MinION Mk1B

1. Connect the MinION to the computer USB Type 3.0 port using the USB Type-A cable provided.
2. Insert the Configuration Test Cell (CTC) as shown.
3. Close the MinION lid.

Note: The MinION lights and fan will now power on after you connect your computer.

4

Perform a **hardware check**

A **hardware check** is required before performing your first MinION Mk1B sequencing run. To run a hardware check, follow the on-screen instructions in MinKNOW, then follow the instructions below. You will require your CTC.

1. In the MinKNOW software, the flow cell status indicator will change colour from grey to white on inserting the MinION CTC into the MinION Mk1B.
2. Click the white box under the MinION positions panel. This will change the colour of the flow cell status indicator on the MinKNOW hardware check panel to dark blue.
3. Press **Start** in the bottom right.
4. Check the flow cell position in MinKNOW shows a  to pass the hardware check.
5. Remove the CTC after you complete the hardware check.

Note: If your hardware check fails, see Support in Additional information.

5


Discover the **Nanopore Community**

community.nanoporetech.com

Ensure the success of your nanopore sequencing project and stay up-to-date with the latest technology and protocol updates.

Follow the link to Getting Started with your experiments to begin sequencing:
community.nanoporetech.com/to/minion-lab

Tip: Learn how to analyse your nanopore data at:
nanoporetech.com/analyse

Additional information

Warranty

A license and warranty can be purchased for your device here:

store.nanoporetech.com/device-warranty.html

Flow cell warranty: community.nanoporetech.com/to/warranty

Recycle used flow cells

Oxford Nanopore is committed to environmental sustainability. You can help by sending your flow cells for recycling.

Find out how: community.nanoporetech.com/support/returns

Place your next order

Buy more consumables at the Oxford Nanopore Store:

store.nanoporetech.com

Documentation

Documentation for your device is available on the Nanopore Community:

community.nanoporetech.com/docs

Support

For all of your customer and technical support needs, visit:

community.nanoporetech.com/support

Technical specification

	MinION Mk1B
Model number	MIN-101B
Supply voltage (V)	5 DC
Maximum rated current (A)	1.0
Maximum rated power (W)	5
Size (H x W x D) (mm)	23 x 33 x 105
Weight (g)	100
Installation ports	1 x USB Type-B
Software installed	MinION driver*
Compute specification	N/A
Environmental conditions	Functional range of electronics is within environmental temperatures of +5°C to +40°C. Use within 30%-75% relative non-condensing humidity limits. Designed to sequence in environmental temperatures of +18°C to +25°C. Intended for indoor use. Can be used up to altitudes of 2,000 m. The device has a Pollution Degree 2.

*Device drivers are used to correctly configure a computer that the device is plugged in to.

Oxford Nanopore Technologies

phone +44 (0)845 034 7900

email support@nanoporetech.com

X @nanopore

www.nanoporetech.com


Oxford Nanopore Technologies, the Wheel icon, MinION, Flongle, and MinKNOW are registered trademarks of Oxford Nanopore Technologies in various countries. © 2024 Oxford Nanopore Technologies. All rights reserved. Oxford Nanopore Technologies products are not intended for use for health assessment or to diagnose, treat, mitigate, cure, or prevent any disease or condition.

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