



Q-Line sequencing software user guide

V Q_SSU_revB_03Feb2025

Contents

Introduction to the Q-Line sequencing software

1. Introduction

Log in and launch the application

2. Log in and launch the application
3. Changing your password

Create a sample sheet

4. Create a sample sheet

Set up and run an assay

5. Set up and run an assay

Run complete

6. Run complete
7. Data location

Technical run report

8. Technical run report

1. Introduction

Introduction to the sequencing software

The sequencing software carries out several core tasks: data acquisition, real-time analysis and feedback, basecalling, data streaming, controlling the device, and ensuring that the platform chemistry is performing correctly to run the samples. The sequencing software takes the raw data and converts it into reads by recognition of the distinctive change in current that occurs when a DNA strand enters and leaves the pore. It then basecalls the reads and writes out the data into FASTQ and BAM files.

Computer specification

Please refer to the **Q system GridION IT requirements** document for details of the IT setup required to run a GridION Q.


2. Log in and launch the application

- 1 **Click on the user account. You will be asked for your password.**

If you have created other user accounts besides GridION user, you will see them in a list.




You will only see local accounts, not Active Directory (AD) accounts. To log in using an AD account, click **“Not listed?”**.



Then, enter the account name without the domain name (e.g., if the domain account name is `John.Smith@example.local`, enter John.Smith in the username prompt).

Note: If you are rebooting your device, you will be asked for your device passphrase, which differs from your user account password. For more information on how to set up your user account and change your device passphrase, see the "User management - local user" of the Q-line software installation and maintenance guide.

- 2 Enter the password. If you have not changed the password for GridION user, the password will be “grid”.**




3 In the sidebar, click the Sequencing Software icon.



Once you open the Sequencing Software, the Overview page will be displayed. The overview


shows the five flow cell positions and their current status.




From this page, you can start running an assay.

There are four buttons on the left side:


- **Overview:** The overview includes the five GridION positions and their current status. From this page, it is possible to set up an assay through the Set up run buttons.
- **Run setup:** The run setup opens to the assay selection. It shows the five stages of setting up a run.



- **Run log:** The run log contains all your completed runs.



- **System:** This section includes information on device storage and language settings. If you have Laboratory Manager or IT Administrator permissions, you will also see Assay Configuration and API Key Management.



From this page, you can begin the following tasks: set up a run, view in-progress or completed runs, manage the device, or view details for live positions.


There are four buttons on the left side:

Overview

The overview includes the five GridION positions, and you can interact on this page to set up your experiment.


Run Setup

Run setup opens up to the assay selection. You can see the five run setups available: Assay selection, flow cell check, sample sheet, prime and load, review and start.




Run log

The run log contains all of the runs you have completed.




System

This section includes: Access device information, Language settings and, if you have Laboratory Manager or IT Administrator permissions, Assay configuration and API key management.




3. Changing your password


- 1 Search "Settings" in the search bar and click on Settings.



- 2 View the GridION username and the current password by clicking the chevron next to the password.



3 Change your password by entering your current and new passwords in the relevant boxes. Make sure you follow the password requirements listed below.





Password requirements

Note: IT Administrators will be able to change these complexity requirements.

- You can use numbers, letters, upper/lower case and punctuation.
- The minimum password length is eight characters.
- You cannot have a password that is a palindrome (reads the same backwards and forwards, e.g. "racecar").
- The new password cannot be the same as the previous one, even if you have changed the cases.
- The new password cannot be the rotated version of the old password (e.g. "billy" and "illyb").
- Multiple character groups must be included in the new password.
- The password must contain at least one character from three of four character groups. These groups are:
 - Digits
 - Lower case alphabetic characters from the current system locale
 - Upper case alphabetic characters from the current system locale
 - Other characters, e.g. punctuation marks and symbols such as * % ! ?

For local users where there is no distinction between upper and lower case (i.e. most Asian scripts), the two groups are treated as one.

You can also click the button next to New Password to create a strong password.



4. Create a sample sheet

When setting up an assay run, you will need to provide a sample sheet detailing the barcodes used, the corresponding sample IDs, and their sample type. This is a vital traceability link, and strict validation rules are applied to minimise the risk of data entry errors.

You can manually create sample sheets on the device using LibreOffice Calc or text editor, available on the GridION. You can also create the sample sheet off-device and transfer it onto the device through drive mounting (to do this, speak to your IT department). You can also transfer your sample sheets via a USB stick (to do this, enable the USB mount as described in the "Managing data" section of the Q-Line Sequencing Software installation and maintenance guide).

Alternatively, you can use a laboratory information management system (LIMS) to add information to your sample sheet automatically. The GridION and LIMS must be configured before library preparation, as detailed in the "Configuring LIMS integrations" section of the Q-Line Sequencing Software installation and maintenance guide.

Manual input

For laboratories without a LIMS, sample sheets can be written manually using Microsoft Excel, LibreOffice Calc, or a text editor. Whichever is used, it is vital to save the file as in the .csv format.

LibreOffice Calc is installed on your GridION. The following steps describe how to use it to generate a sample sheet.

- 1 **Open the template for the assay sheet in LibreOffice Calc.**
- 2 **Edit the template for the samples you are running.**
- 3 **Save the sample sheet as a .csv file (not ODS template) in your local folder.**
- 4 **Import the sample sheet when setting up an assay, as described in the "Set up and run and assay" section of this guide.**

Template format: available as templates (.ods files) in your GridION

| Keys | Description |
|---------------------------|---|
| v1 (above header section) | Version of the sample sheet format (Do not change this) |
| assay | Assay key installed on the UI and its version. Must be in the format key:version where key and version are defined in the assay definition file |
| library_id | Provided by you to ID the library for sequencing (Any of <code>a-z A-Z 0-9 \.\-_]+[,\s]*\$)</code> Important: Do not add any personally-identifiable information in this field. |
| sample_count | Number of samples used. This must match the number of barcodes used, including samples and controls. Note: Using the template will automatically calculate this from the barcode/sample section. |
| created_by | (Optional) User who prepared the library (Any of <code>a-z A-Z 0-9 \.\-_]+[,\s]*\$)</code> |
| lot | (Optional) The kind of item, and the unique ID for lot used. (In separate cells on the same row, see image below for an example) |

Barcode/sample section

barcode

The barcode used during library preparation for the sample. **Important:** If you are not using barcoding, you will not have the barcode column and will only have a sample count of 1.

sample_id

Provided by you to ID individual input samples. It must be unique within the assay. **Important:** Do not add any personally-identifiable information in this field.

sample_type

The type of sample. For more information on acceptable fields, see the "Configuring an assay definition file for your own assay" section of the Q-Line Sequencing Software Installation and maintenance guide. This field will depend on the assay run and its definitions, such as sample and control sample.

Examples

The table below shows the sample sheet for an assay using the Ligation Sequencing Kit without barcoding. Important note: empty rows below `v1` and `sample_count` rows must be conserved in your own sample sheets. A failure to do so will invalidate the document.

| v1 | |
|--------------|----------------|
| assay | ligation-kit09 |
| library_id | My_Library_ID |
| sample_count | 1 |
| sample_type | sample_id |
| sample | Test_01 |

Below is an example sample sheet for the Rapid Barcoding Kit:


| | | |
|--------------|-----------------------------|------------|
| v1 | | |
| assay | rapid-barcode-96-kit10:v1.0 | |
| library_id | My_Library_ID | |
| sample_count | 5 | |
| barcode | sample_type | sample_id |
| barcode01 | no_template_control | Control_01 |
| barcode02 | no_template_control | Control_02 |
| barcode03 | sample | Test_01 |
| barcode04 | sample | Test_02 |
| barcode05 | sample | Test_03 |

5. Set up and run an assay

This section explains how to set up and conduct your assay. We will use the Rapid Barcoding Experiment 01-96 assay as an example, and will guide you through selecting your assay, performing a flow cell check, adding a sample sheet, priming and loading your flow cell, and initiating your run.


1 Click **"Overview"** to see the five flow cell positions. Then click **"Set up run"**.

You can also click **Run setup** in the sidebar. **Note:** The **Set up run** option is only available if flow cells are inserted and available for use.




2 Choose the Rapid Barcoding Experiment 01-96 assay, and click Continue.


You can also add custom assays to this list; the method to generate a custom assay is described in the "Configuring an assay definition file for your own assay" of the Q-Line Sequencing Software Installation and maintenance guide.



3 Click Check flow cell. This will initiate the check for whether the flow cell is usable for the assay.




A status message will appear under the flow cell position: "Checking flow cell...". You can also cancel it by clicking **Cancel**.



4 Once the flow cell check has finished, the software will report whether the flow cell has sufficient health to run the assay. Click Continue to move to the next step.

A failed flow cell check will prevent the start of a run. However, it is possible to override a low flow cell health check with certain assays, including most default assays and custom assays if configured that way. If your assay allows the overriding of low flow cell health, the number of available pores is displayed to enable you to make an informed decision.

In some cases, the low flow cell health might have been caused by poor thermal contact between the position and the flow cell. In these cases, on-screen instructions will direct you to reinsert the flow cell and try again.





5 Select the sample sheet that you want to use. You can only select sample sheets that have a Validation status of Pass. If the sample sheet you want to use has a different status, click View sheet for more information. You may need to use a different flow cell.

Tip: If you have a barcode scanner plugged into the GridION, you can scan a barcode for the Library ID of the sample sheet you want to use and filter the list to only display the matching sample sheet.

If your sample sheet is not already displayed, you can create and upload one. The method for creating a sample sheet can be found in the "Create a sample sheet" section of this guide.


To manually upload your sample sheet, click **Import**.

The sequencing software will check whether the sample sheet meets the requirements. If no issues have been found, click **Add sample sheet**.




If your sample sheet has an issue, an error message will pop up, advising you which part of the sample sheet does not match the required criteria.

Occasionally, when adding the sample sheet, the progress bar may not move (as shown in the image below). However, clicking Add sample sheet will enable the sample sheet to be added.



- 6 When the sample sheet has been imported, a green tick will show it has passed validation. You can also check your sample sheet by clicking the View sheet button.



- 7 When you are happy with your sample sheet selection, click Continue.

- 8 Prime and load your flow cell. You can see detailed instructions for this in the View full instructions tab.

The screenshot shows a software interface for a flow cell setup. On the left is a sidebar with icons for Overview, Run setup, Run log, and System. The main window has a progress bar with five steps: 1 Assay selection, 2 Flow cell check, 3 Sample sheet, 4 Prime and load (active), and 5 Review and start. The 'Prime and load' section contains an 'ACTION REQUIRED' box with instructions: 'Prime flow cell in position X1', 'Load prepared library into flow cell', and a 'View full instructions' button. Below this is a grid of five flow cell positions (X1-X5) with their respective statuses: X1 (Prime flow cell and load library), X2 (Flow cell incompatible), X3 (In use), X4 (Flow cell incompatible), and X5 (Flow cell detected). At the bottom are 'Cancel', 'Back', and 'Continue' buttons.

Run setup grid

1 Assay selection 2 Flow cell check 3 Sample sheet 4 Prime and load 5 Review and start

Prime and load

ACTION REQUIRED

Prime flow cell in position X1
It is recommended to prime and load on the device.

Load prepared library into flow cell
Add 75 µl of sample library.
Click 'Continue' when done.


[View full instructions](#)

| X1 | X2 | X3 | X4 | X5 |
|----------------------------------|------------------------|--------|------------------------|--------------------|
| Prime flow cell and load library | Flow cell incompatible | In use | Flow cell incompatible | Flow cell detected |


[Cancel](#) [Back](#) [Continue](#)

System


- 9 Once you have loaded your sample, click Continue.




10 You can now review the assay, libraryID, number of samples, and time to report and check your sample sheet.



- 11 If you are happy with the settings, click Run assay. A small prompt will appear in the bottom right corner stating that the run has started. If you did not mean to start it, click Undo within 30 seconds, and the run will stop.



- 12 Click the View run button to see the status of your run.



The **View run** button will open the Run log screen for a particular position. There, four panel are displayed: Run status, Run summary, Run reports and Samples.

- **Run status:** displays the Initial health check and the Initial sequencing performance. Both reports on potential issues with the run:
 - **Initial health check:** reports on issues introduced between the flow cell check and the sequencing.
 - **Initial sequencing performance:** reports on potential issue with the library.
- **Run summary:** recapitulates the configuration of the run.
- **Run reports:** provides a list of reports available.
- **Samples:** If you have barcodes, you can also check whether individual samples barcodes are found.

Note that this information is also available in the Overview page.



Overview

Run setup

Run log

System

Run log / Run for Rapid Barcoding Experiment 01-96 23-08-24 16:38

grid

RUN FOR LIBRARY ID

RBK-LIBRARY-ID In progress

Cancel run

Run status



TIME TO RUN COMPLETION

About 24 hours

RUN HEALTH

Initial health check No issues found

Initial sequencing performance No issues found

Run summary

| | |
|--------------|----------------------------------|
| Library ID | rbk-library-id |
| Assay type | Rapid Barcoding Experiment 01-96 |
| Sample count | 8 |
| Position | X1 |
| Operator | grid |
| Start time | 23-08-24 16:38 |

Run reports

| Report | Status | Progress |
|--------|--------|----------|
|--------|--------|----------|

| | | |
|----------------------|--|----------------------|
| Technical run report | | View |
|----------------------|--|----------------------|

Samples


[View sample sheet](#)

Samples found

2 / 8

6. Run complete

After the run finishes, view the run report by clicking Run log. It provides details on the run summary, run status, run reports, and sample information.



Run log / Run for Rapid Barcoding Experiment 01-96 23-08-24 16:38 grid

RUN FOR LIBRARY ID
RBK-LIBRARY-ID Completed

Run status

TIME TO RUN COMPLETION
Finished

RUN HEALTH

- Initial health check No issues found
- Initial sequencing performance No issues found

Run summary

| | |
|--------------|----------------------------------|
| Library ID | rbk-library-id |
| Assay type | Rapid Barcoding Experiment 01-96 |
| Sample count | 8 |
| Position | X1 |
| Operator | grid |
| Start time | 23-08-24 16:38 |

Run reports

| Report | Status | Progress |
|----------------------|--------|----------------------|
| Technical run report | | View |

Samples [View sample sheet](#)

Samples found
4 / 8

| Sample ID | Type | Sample status | Sample reports |
|-----------------|---------------------|---------------|----------------|
| sample1 BR20 | sample | Not found | Not available |
| sample2 BR21 | sample | Found | Not available |
| sample3 BR22 | no_template_control | Not found | |
| sample4 BR23 | sample | Found | Not available |
| sample5 BR24 | sample | Not found | Not available |
| sample6 BR25 | sample | Found | Not available |
| sample7 BR26 | sample | Found | Not available |
| sample8 BR27 | sample | Not found | Not available |

1 Click View under the run reports section to access the technical run report.

The report contains information relating to the sequencing performance of your run, including a variety of graphs and statistics. These are useful for troubleshooting and quality control.



Overview

Run setup

Run log

System

RUN FOR LIBRARY ID

RBK-LIBRARY-ID Completed

Run status



TIME TO RUN COMPLETION

Finished

RUN HEALTH

Initial health check No issues found

Initial sequencing performance No issues found

Run summary

| | |
|--------------|----------------------------------|
| Library ID | rbk-library-id |
| Assay type | Rapid Barcoding Experiment 01-96 |
| Sample count | 8 |
| Position | X1 |
| Operator | grid |
| Start time | 23-08-24 16:38 |

Run reports

| Report | Status | Progress | |
|----------------------|--------|----------|----------------------|
| Technical run report | | | View |

Samples

[View sample sheet](#)

Samples found

4 / 8

| Sample ID | Type | Sample status | Sample reports |
|-----------------|---------------------|------------------------|----------------|
| sample1 BR20 | sample | Not found | Not available |
| sample2 BR21 | sample | Found | Not available |
| sample3 BR22 | no_template_control | Not found | |
| sample4 BR23 | sample | Found | Not available |
| sample5 BR24 | sample | Not found | Not available |
| sample6 BR25 | sample | Found | Not available |
| sample7 BR26 | sample | Found | Not available |
| sample8 BR27 | sample | Not found | Not available |

You can save the final technical run report as a HTML file by clicking the **Save a copy** button in the top right corner of the screen.

Technical run report Save a copy Close

GridION Mk1 (INT-GRD-1644) Final report

23 Aug 24, 16:38 — 24 Aug 24, 16:41 · no_group · no_sample · X1
Protocol run ID: eba9736b-4312-4732-a0f5-0556d2aaa13d

[Run summary](#) | [Run configuration](#) | [Sequence output](#) | [Run health](#) | [Run log](#)

Run summary

DATA OUTPUT

| | | | |
|-----------------|---------|---------------|----------|
| Estimated bases | 4.05 Gb | Data produced | 13.21 GB |
| Reads generated | 11.25 M | Estimated N50 | 334 |

RUN DURATION

Run time: 24 hrs 0 mins / 24 hrs 0 mins

Elapsed time: Run limit:

Run status: **FINISHED** · Target runtime has been reached

[View unit abbreviations used in this report](#)

BASECALLING

| | | | |
|--------------|------|-------------------------------|---------------|
| Reads called | 100% | Bases called (min Q score: 9) | |
| | | 3.66 Gb (Pass) | 5.6 Mb (Fail) |

Run configuration

RUN SETUP

| | |
|----------------|---------------|
| Flow cell type | FLO-MIN106 |
| Flow cell ID | FAW97034 |
| Kit type | SQK-RBK110-96 |

RUN SETTINGS

| | |
|--------------------------|--------|
| Run limit | 24 hrs |
| Active channel selection | On |

DATA OUTPUT SETTINGS

| | |
|----------------------|--|
| FAST5 output | Off |
| FASTQ output | gzip_compress |
| FASTQ reads per file | 4000 |
| BAM output | Off |
| Bulk file output | Off |
| Data location | /data/output/sequences/./no_group/no_sample/20240823_1638_X1_FAW97034_96 |

Note: Additional reports (in PDF or HTML format) will figure in the Sample ID list if an analysis was available for your assay.

Example: Instrument qualification run report

Run log / Run for Instrument Qualification 04-07-24 14:08 grid

RUN FOR LIBRARY ID
IQ-LIBRARY-ID Completed

Run status

✓

TIME TO RUN COMPLETION
Finished

RUN HEALTH

Initial health checkNo issues found

Initial sequencing performanceNo issues found

Run summary

| | |
|--------------|--------------------------|
| Library ID | iq-library-id |
| Assay type | Instrument Qualification |
| Sample count | 8 |
| Position | X4 |
| Operator | grid |
| Start time | 04-07-24 14:08 |

Run reports

| Report | Status | Progress | |
|----------------------|--------------------|---|----------------------|
| Workflow report | ● Report available | <div style="width: 100%; height: 10px; background-color: #2c3e50;"></div> | View |
| Workflow report pdf | ● Report available | <div style="width: 100%; height: 10px; background-color: #2c3e50;"></div> | PDF |
| Technical run report | | | View |

Samples

Samples found

Samples failed

[View sample sheet](#)

Workflow report pdf Save a copy X Close

1 of 2 Page Width

Instrument Qualification

Oxford Nanopore Technologies/Diagnostics Device Software System Test

REPORT GENERATED ON
2024-07-04 14:23:30.762303

ANALYSIS VERSION
v1.0.0

DEVICE
GridION

STARTED AT
2024-07-04T14:08:14.264650+00:00

STARTED BY
grid

DEVICE SERIAL NUMBER
INT-GRD-1644

DEVICE TYPE
GridION Mk1

LIBRARY ID
iq-library-id

SEQUENCING SOFTWARE VERSION
23.06.19

FLOW CELL ID
FAW99036

WORKFLOW BROKER
5.1.13

Instrument qualification status ❌ Fail

7. Data location

If data offload is configured, sequencing and analysis data and reports are moved to:

```
/data/offload/<assay_run_id>/
```

From here, they are moved to the network drive location set by the IT administrator. For data offload instructions, refer to the “Managing data” section of the Q-Line sequencing software installation and maintenance guide.

Sequencing data and reports are located in:

```
/data/offload/<assay_run_id>/sequencing
```

The file structure within the sequencing folder is as follows:

```
|-- /data/offload/<assay_run_id>/sequencing
|   |-- barcode_alignment.tsv
|   |-- final_summary.txt
|   |-- pore_activity.csv
|   |-- report.html
|   |-- report.json
|   |-- report.md
|   |-- sample_sheet.csv
|   |-- sequencing_summary.txt
|   |-- throughput.csv
|   |-- fastq_fail
|   |   |-- barcode01
|   |   |   |-- <flow_cell_id>_fail_barcode01_<id>.fastq.gz
|   |   |-- barcode02
|   |   |   |-- <flow_cell_id>_fail_barcode02_<id>.fastq.gz
|   |   |-- unclassified
|   |   |   |-- <flow_cell_id>_fail_unclassified_<id>.fastq.gz
|   |-- fastq_pass
|   |   |-- barcode01
|   |   |   |-- <flow_cell_id>_pass_barcode01_<id>.fastq.gz
|   |   |-- barcode02
|   |   |   |-- <flow_cell_id>_pass_barcode02_<id>.fastq.gz
|   |   |-- unclassified
|   |   |   |-- <flow_cell_id>_pass_unclassified_<id>.fastq.gz
|   |-- other_reports
|   |   |-- pore_scan_data.csv
```

The output location for analysis data and reports is:

```
/data/offload/<assay_run_id>/analysis
```

The file structure within the analysis folder is as follows. There is a set of output files for every sample/replicate run in the assay.

```
|-- /data/offload/<assay_run_id>/analysis
|   |-- <output_file>.vcf
|   |-- <output_file>.bam
|   |-- <output_file>.bam.bai
|   |-- <output_file>.html
|   |-- <output_file>.pdf
|   |-- execution
|   |   |-- report.html
|   |   |-- timeline.html
|   |   |-- trace.txt
```

8. Technical run report

Technical run report overview

Run reports contain information about the sequencing run and include performance graphs. These graphs are automatically generated at the end of a run or when you click the **View** button in the **Run report** panel of the **Run log** tab.

The screenshot displays the 'Run log' interface for a completed run. The main content area is titled 'RUN FOR LIBRARY ID RBK-LIBRARY-ID' with a 'Completed' status indicator. The interface is divided into several sections:

- Run status:** Shows 'TIME TO RUN COMPLETION' as 'Finished' with a progress bar. Below this, 'RUN HEALTH' is displayed with two green bars: 'Initial health check' (No issues found) and 'Initial sequencing performance' (No issues found).
- Run summary:** A table listing key run parameters:

| | |
|--------------|----------------------------------|
| Library ID | rbk-library-id |
| Assay type | Rapid Barcoding Experiment 01-96 |
| Sample count | 8 |
| Position | X1 |
| Operator | grid |
| Start time | 23-08-24 16:38 |
- Run reports:** A table with columns 'Report', 'Status', and 'Progress'. A 'Technical run report' is listed, with a 'View' button highlighted by a red box.
- Samples:** A section titled 'Samples found 4 / 8' with a 'View sample sheet' link. Below is a table of sample results:

| Sample ID | Type | Sample status | Sample reports |
|-----------------|---------------------|---------------|----------------|
| sample1 BR20 | sample | Not found | Not available |
| sample2 BR21 | sample | Found | Not available |
| sample3 BR22 | no_template_control | Not found | Not available |
| sample4 BR23 | sample | Found | Not available |
| sample5 BR24 | sample | Not found | Not available |
| sample6 BR25 | sample | Found | Not available |
| sample7 BR26 | sample | Found | Not available |
| sample8 BR27 | sample | Not found | Not available |

The run report includes panels for the Run summary, Run configuration, Sequence output, Alignment (only if live alignment is part of your assay), Run health and Run log. It is interactive: sections of interest can be expanded and navigated while troubleshooting suggestions are made available for performance enhancement.

Run summary

This is an overview of the sequencing run, including output and basecalling results and the duration of the run.

^ Run summary

DATA OUTPUT

Estimated bases

8.58 Gb

Data produced

10.64 GB

Reads generated

3.39 M

Estimated N50

5.78 kb

BASECALLING

Reads called

100%

Bases called (min Q score: 9)

6.95 Gb

Pass

1.06 Gb

Fail

RUN DURATION

Run time

24 hrs 0 mins / 24 hrs 0 mins

Elapsed time Run limit

Run status

FINISHED - Target runtime has been reached

[View unit abbreviations used in this report](#)

Run configuration

A detailed overview of the settings selected for the sequencing run and data output, including the software versions used to sequence the data. It also includes further information about the modified bases and basecalling models used, expansion kits used and read-splitting preferences.

^ Run configuration

^ RUN SETUP

| | |
|----------------|---------------|
| Flow cell type | FLO-MIN106 |
| Flow cell ID | FAZ39326 |
| Kit type | SQK-RBK110-96 |

^ RUN SETTINGS

| | |
|--|---|
| Run limit | 24h |
| Active channel selection | On |
| Pore scan freq. | 1.5 hrs |
| Reserved pores | On |
| Minimum read length | 200 bp |
| Read splitting | On |
| Basecalling | High-accuracy model, 450 bps |
| Modified basecalling | Off |
| Trim barcodes | Off |
| Mid-read barcode filtering | On |
| Overridden value of minimum mid-read barcoding score | 60 |
| Alignment reference | /data/references/lambda_75_7a991a.fasta |

^ DATA OUTPUT SETTINGS

| | |
|----------------------|---|
| FAST5 output | Off |
| FASTQ output | gzip_compress |
| FASTQ reads per file | 4000 |
| BAM output | Off |
| Bulk file output | Off |
| Data location | /data/output/sequences/./no_group/no_sample/20240618_1446_X3_FAZ39326_803_990b5 |


^ SOFTWARE VERSIONS

| | |
|---------------|----------|
| MinKNOW | 23.06.16 |
| Bream | 7.6.6 |
| Configuration | 5.6.5 |
| Guppy | 6.6.3 |
| MinKNOW Core | 5.6.3 |

Sequence output

This section includes a more detailed view of the sequencing output, such as read lengths, barcodes detected and quality score. This section of the Technical run report displays Read lengths, Barcodes, Cumulative output and Quality scores data.

Read lengths This graph shows the total number of bases against the read length. Any outliers are displayed in the table below, displaying read length and aggregated reads.



Barcoded reads

This section will be visible only if barcoding is enabled in the assay. It will then show the total number of bases for each barcode.

^ BARCODED READS

The total number of bases for each barcode is calculated and displayed below.

| Total bases (Gb) | Barcodes |
|-----------------------------|-----------------------------|
| 0-1 | ONQ-test-001 (Reads: 23549) |
| | ONQ-test-002 (Reads: 35496) |
| | ONQ-test-003 (Reads: 20496) |
| | ONQ-test-004 (Reads: 27890) |
| | ONQ-test-005 (Reads: 33134) |
| | ONQ-test-006 (Reads: 20655) |
| | ONQ-test-007 (Reads: 18822) |
| | ONQ-test-008 (Reads: 17256) |
| | ONQ-test-009 (Reads: 14849) |
| | ONQ-test-010 (Reads: 17284) |
| | ONQ-test-011 (Reads: 31763) |
| | ONQ-test-012 (Reads: 25841) |
| | ONQ-test-013 (Reads: 25695) |
| | ONQ-test-014 (Reads: 16672) |
| | ONQ-test-015 (Reads: 12485) |
| | ONQ-test-016 (Reads: 24885) |
| | ONQ-test-017 (Reads: 18270) |
| | ONQ-test-018 (Reads: 26146) |
| | ONQ-test-019 (Reads: 25975) |
| | ONQ-test-020 (Reads: 27983) |
| | ONQ-test-021 (Reads: 25919) |
| | ONQ-test-022 (Reads: 17141) |
| | ONQ-test-023 (Reads: 15591) |
| | ONQ-test-024 (Reads: 16890) |
| | ONQ-test-025 (Reads: 30629) |
| | ONQ-test-026 (Reads: 20670) |
| | ONQ-test-027 (Reads: 29451) |
| | ONQ-test-028 (Reads: 37344) |
| | ONQ-test-029 (Reads: 35602) |
| ONQ-test-030 (Reads: 30579) | |
| ONQ-test-031 (Reads: 24205) | |
| ONQ-test-032 (Reads: 20860) | |
| ONQ-test-033 (Reads: 25382) | |
| ONQ-test-034 (Reads: 20483) | |
| ONQ-test-035 (Reads: 21515) | |
| ONQ-test-036 (Reads: 17710) | |
| ONQ-test-037 (Reads: 35562) | |
| ONQ-test-038 (Reads: 24063) | |
| ONQ-test-039 (Reads: 26001) | |
| ONQ-test-040 (Reads: 20218) | |
| ONQ-test-041 (Reads: 25747) | |
| ONQ-test-042 (Reads: 20724) | |
| ONQ-test-043 (Reads: 23405) | |
| ONQ-test-044 (Reads: 27288) | |
| ONQ-test-045 (Reads: 34220) | |
| ONQ-test-046 (Reads: 27375) | |
| ONQ-test-047 (Reads: 21834) | |
| ONQ-test-048 (Reads: 16034) | |
| ONQ-test-049 (Reads: 26690) | |
| ONQ-test-050 (Reads: 33880) | |
| ONQ-test-051 (Reads: 25087) | |
| ONQ-test-052 (Reads: 32458) | |
| ONQ-test-053 (Reads: 27059) | |
| ONQ-test-054 (Reads: 14978) | |
| ONQ-test-055 (Reads: 48230) | |
| ONQ-test-056 (Reads: 36834) | |
| ONQ-test-057 (Reads: 39493) | |
| ONQ-test-058 (Reads: 12140) | |
| ONQ-test-059 (Reads: 32303) | |
| ONQ-test-060 (Reads: 25732) | |
| ONQ-test-061 (Reads: 18747) | |
| ONQ-test-062 (Reads: 17926) | |
| ONQ-test-063 (Reads: 31772) | |
| ONQ-test-064 (Reads: 27197) | |
| ONQ-test-065 (Reads: 35834) | |
| ONQ-test-066 (Reads: 28975) | |
| ONQ-test-067 (Reads: 19185) | |
| ONQ-test-068 (Reads: 26458) | |
| ONQ-test-069 (Reads: 25939) | |
| ONQ-test-070 (Reads: 24657) | |
| ONQ-test-071 (Reads: 26433) | |
| ONQ-test-072 (Reads: 33964) | |
| ONQ-test-073 (Reads: 36936) | |
| ONQ-test-074 (Reads: 40260) | |
| ONQ-test-075 (Reads: 15649) | |
| ONQ-test-076 (Reads: 18256) | |
| ONQ-test-077 (Reads: 31932) | |
| ONQ-test-078 (Reads: 19693) | |
| ONQ-test-079 (Reads: 30818) | |
| ONQ-test-080 (Reads: 29165) | |
| ONQ-test-081 (Reads: 39965) | |
| ONQ-test-082 (Reads: 24517) | |
| ONQ-test-083 (Reads: 26791) | |
| ONQ-test-084 (Reads: 18037) | |
| ONQ-test-085 (Reads: 24352) | |
| ONQ-test-086 (Reads: 22109) | |
| ONQ-test-087 (Reads: 44907) | |
| ONQ-test-088 (Reads: 53720) | |
| ONQ-test-089 (Reads: 32214) | |
| ONQ-test-090 (Reads: 34864) | |
| ONQ-test-091 (Reads: 19972) | |
| ONQ-test-092 (Reads: 42084) | |
| ONQ-test-093 (Reads: 23724) | |
| ONQ-test-094 (Reads: 24443) | |
| ONQ-test-095 (Reads: 26309) | |
| ONQ-test-096 (Reads: 25199) | |

Cumulative output

The graphs show the total output of the bases and reads sequenced during the experiment.


CUMULATIVE OUTPUT

The cumulative output shows the total amount of bases or reads sequenced over time by your device.

Bases

Legend


- Estimated** (Blue line): Predicted total number of bases, prior to basecalling
- Passed** (Green line): Bases equal to or above the quality score threshold.
- Failed** (Red line): Bases below the quality score threshold.



Reads

Legend

- Total** (Blue line): Total number of reads, including passed, failed and skipped.
- Passed** (Green line): Reads equal to or above the quality score threshold.
- Failed** (Red line): Reads below the quality score threshold.
- Skipped** (Grey line): Reads that will not be basecalled. Post run basecalling is possible.



Quality score


The quality score is calculated as basecalling is performed on your device. Reads that fall below the minimum value of 9 will be classified as failed reads. You can alter the accepted minimum quality score in the assay configuration file.

^ QUALITY SCORE

The quality score is calculated as basecalling is performed on your device. Reads that fall below the minimum value of 9 will be classified as failed reads. You can alter the accepted minimum quality score in MinkNOW.

Legend

- Mode**
The most frequent quality score of reads in the run.
- ▭ Spread**
The spread of quality scores, found by calculating full width half maximum.
- - - Min. quality score**
Minimum quality score to be accepted as a passed read.



^ Troubleshooting

Quality score low

This can be due to the translocation speed being out of the accepted range, which can correlate to low quality scores. If you see that the translocation speed is out of the accepted range in the below graph, please see the Flow Cell refuelling page linked [here](#) for further troubleshooting.

Alignment

This section will figure in your Technical run report if an alignment file was provided and Alignment enabled in the assay configuration file.

Reference alignment: This section will display the total number of passed reads aligned to each uploaded reference target, calculated and displayed.

Bed regions: If there is corresponding .bed region information for aligned references, it will be shown here.

^ Alignment

^ Reference Alignment

The total number of passed reads aligned to each uploaded reference target is calculated and displayed below.

| Total reads | Reference target |
|-------------|------------------|
| 0–1 k | <i>None</i> |
| 1 k–10 k | <i>None</i> |
| 10 k–50 k | <i>None</i> |
| 50 k–100 k | <i>None</i> |
| 100 k–500 k | <i>None</i> |
| 500 k+ | <i>lambda</i> |

^ Bed Regions

Corresponding bed region information for aligned references

| Coordinate | Name | Read Count |
|------------|------|------------|
|------------|------|------------|

Run health


This section provides a detailed view of the Pore activity and Pore scan results throughout the run, along with graphs displaying the translocation speed and temperature.

PORE ACTIVITY

The Pore activity graph shows the performance of your sample as it is being sequenced during a run.

Legend

- Sequencing
Pore currently sequencing
- Adapter
Pore currently sequencing adapter
- Pore available
Pore available for sequencing
- Unavailable
Pore unavailable for sequencing
- Active feedback
Channel ejecting analyte
- No pore
No pore detected in channel
- Out of range-high
Current is positive but unavailable for sequencing
- Out of range-low
Current is negative but unavailable for sequencing
- Multiple
Multiple pores detected. Unavailable for sequencing.
- Saturated
The channel has switched off as current levels exceed hardware limitations
- Zero
Pore currently unavailable for sequencing
- Channel disabled
Channel is disabled and awaiting another pore scan
- Unclassified
Pore status unknown



Troubleshooting

General


Some commonly seen issues are excess pores classified as Recovering, Open Pore, or Free Adapter. To find out what advice is applicable for your run, visit the [user guide](#).

PORE SCAN

A Pore scan is performed at configurable time intervals to determine the current status of pores within channels on a Flow Cell. For this run a Pore scan is performed every 1.5 hrs.

Legend

- Pore available
Pore in channel available for sequencing
- Reserved pore
Pore in reserve, will return to available when required
- Unavailable
Pore inhibited from sequencing
- Saturated
Possible contamination in the sample
- Zero
No current is passing through this pore, possibly due to bubbles on the membrane
- Inactive
Pore no longer suitable for further sequencing



Troubleshooting

High proportion Unavailable
Possible contaminants in library blocking the pore. Consider using the Flow Cell Wash Kit, and reloading a library.

High proportion Inactive
If localised to one area of the Flow Cell, this could indicate that an air bubble has been introduced during the flushing/loading steps. If inactivity is spread across the Flow Cell this could be caused by improper loading of the library, please refer to the [user guide](#) for further support.


Pore scan

TRANSLOCATION SPEED

The translocation speed is the rate at which DNA/RNA travels through pores as it is being sequenced.

Legend

— Median - - - 75% quartile - - - 25% quartile □ Accepted range



Troubleshooting

Low speed

Check that the Flow Cell is within the target temperature range.

Note


Low-quality and short reads are not included in this graph.

TEMPERATURE

The temperature of the Flow Cell over the run time.

Legend

— Measured — Target



Troubleshooting

Out of range

Check that the Flow Cell is correctly seated and firmly pushed down into the device. Ensure ambient temperature is always within the specified range for your device in the [user guide](#).

Air flow should be good but not excessive. Excessive amounts of cool air blowing on the device could prevent it from reaching target temperature.

Run log

The run log includes system messages sent during the sequencing run regarding errors, warnings, and any events.

^ Run log

SYSTEM MESSAGES

System messages are a record of the events that occurred in the time covered by this report.

^ Errors

None

^ Warnings

None

^ Events

Disk space · 18 Jun 24, 14:46

Disk /data has 6630 GB space remaining

Waiting for temperature · 18 Jun 24, 14:46

Waiting up to 300 seconds for temperature to stabilise at 34.0°C

Starting · 18 Jun 24, 14:51

Starting sequencing procedure

Pore scan starting · 18 Jun 24, 14:51

Performing Pore Scan

Pore scan result · 18 Jun 24, 14:53

Pore scan for flow cell FAZ39326 has found a total of 1487 pores. 494 pores available for immediate sequencing

Pore scan starting · 18 Jun 24, 16:24

Performing Pore Scan

Pore scan result · 18 Jun 24, 16:26

Pore scan for flow cell FAZ39326 has found a total of 1441 pores. 489 pores available for immediate sequencing

Pore scan starting · 18 Jun 24, 17:57

Performing Pore Scan