PromethION 2 Integrated

Quick start guide

Pre-install

This Quick start guide sets up your PromethION™ 2 Integrated for use.


See the following:

PromethION 2 Integrated user manual
community.nanoporetech.com/to/p2i

Safety and regulatory information
community.nanoporetech.com/to/safety

For troubleshooting, see the user manual.

PromethION 2 Integrated

What's in the box:

PromethION 2 Integrated

Power cables*

2 x Configuration Test Cells (CTCs)


What you will need:

DisplayPort or HDMI-compatible monitor

Ethernet cable

Wired USB-A keyboard


Wired USB-A mouse

*The PromethION 2 Integrated ships with 5 x power cables (1 US, 1 UK, 1 EU, 1 CN, 1 AUS) for international use.

1

Set up your device

1. **Unpack** your PromethION 2 Integrated device*.
2. Attach the **cables** and **peripherals** in the order shown opposite (a-c).
Note: You are advised to use one external monitor connected to your HDMI or your DisplayPort (a).
3. Connect your power cable to your **power socket** (d) and **power supply** and switch on the mains, as needed.
4. Press the **power switch** (e).

*Ensure 30 cm clearance to front, rear and sides; keep away from direct heat and sunlight. See user manual for detailed advice.

2

Log in

1. **Log in** to your PromethION 2 Integrated
Username: **prom**
Password: **prom**
2. Update your PromethION 2 Integrated **password**.
Your new password must:
 - Have a minimum 8 characters
 - Contain a capital letter and a lower-case letter
 - Contain a number or special character
 - Not be a palindrome (e.g. "radar")
 - Differ from a previous password by more than 2 characters.
3. Click  on the desktop to load MinKNOW.
4. Log in to MinKNOW using your Oxford Nanopore account details.
5. Familiarise yourself with the tutorials (top right) or click **Skip all tutorials**.

3

Update MinKNOW

1. Click **Host settings** in the side bar.
2. Click **Software** in the side bar.
3. Click **Install update**.
4. Select **MinKNOW only** in the pop-up.

Once downloaded and installed, you must now shut down your PromethION 2 Integrated using the instructions in step 4.

4

Power off


1. Click **Host settings**.
2. Click **Shut down** in the main panel.
3. Click **Shut down** in the pop-up.
4. Press the **power switch** and repeat **Log in** from the previous page (step 2).

After repeating step 2, proceed to '**Perform a Hardware check**' (step 5).

5

Perform a Hardware check

1. Lift the lid of the sequencing positions and insert the CTCs into the device as shown here. Close the lid of the device.
2. Click **Start** (side bar).
3. Click **Hardware check** (under Utilities).
4. Click **Select all positions**.
5. Click **Start** (bottom right).
6. Confirm **Selected position(s) passed all checks** is displayed*.
7. Remove CTCs by sliding them towards you and out of the flow cell sockets. Now close the lid of the sequencing positions.

Your PromethION 2 Integrated is now ready to start sequencing!

*If any hardware check fails: 1. Check the CTCs are fully inserted. 2. Repeat the Hardware check. See **Support** in **Additional information** if the problems persist.

6

The Nanopore Community

Discover the Nanopore Community to find protocols, forums, support and downloads for your experiments.
community.nanoporetech.com

TIP: Learn how to analyse your nanopore data at:
nanoporetech.com/analyse

Additional information

- ✓ **Warranty**
A license and warranty can be purchased for your device here:
store.nanoporetech.com/device-warranty.html
Flow cell warranty, visit: community.nanoporetech.com/to/warranty
- ↻ **Recycle used flow cells**
Oxford Nanopore is committed to environmental sustainability. You can help by sending your flow cells for recycling. Find out how: community.nanoporetech.com/support/returns
- 🛒 **Place your next order**
Buy more consumables at the Oxford Nanopore Store:
store.nanoporetech.com
- 📖 **Documentation**
Documentation for your device is available on the Nanopore Community:
community.nanoporetech.com/docs
- 👍 **Support**
For all of your customer and technical support needs, visit:
community.nanoporetech.com/support

Technical specification

Model number	PRO-INT002	
Supply voltage (V)	100-240 AC ± 10% (50-60 Hz)	
Maximum rated current (A)	10	
Maximum rated power (W)	750	
Size (H x W x D) (mm)	180 x 225 x 430	
Weight (kg)	10.6	
Installation ports	1 x Ethernet port 4 x USB-A 3.0 ports	1 x HDMI port 1 x DisplayPort
Software installed	Ubuntu, PromethION OS, MinkNOW	
Compute specification	5.5" Touchscreen 1 x Intel CPU (12-core) 64 GB Memory	15 TB SSD storage 1 x NVIDIA Ampere-series GPU
Environmental conditions	Functional range of electronics is within environmental temperatures of +5°C to +40°C. Users should allow 30 cm clearance to the rear and sides of the device. Designed to sequence in environmental temperatures of +18°C to +25°C. Use within 30%–75% relative non-condensing humidity limits. Intended for indoor use. Can be used up to altitudes of 2,000 m. The device has a Pollution Degree 2.	

Oxford Nanopore Technologies

phone +44 (0)845 034 7900

email support@nanoporetech.com

 oxford-nanopore-technologies

 @nanopore

 @nanoporetech.com


Gosling Building,
Edmund Halley Road,
Oxford Science Park,
Oxfordshire,
OX4 4DQ, UK

www.nanoporetech.com

Oxford Nanopore Technologies, the Wheel icon, MinKNOW and PromethION, are registered trademarks or the subject of trademark applications of Oxford Nanopore Technologies plc in various countries. Information contained herein may be protected by patents or patents pending of Oxford Nanopore Technologies plc. All other brands and names contained are the property of their respective owners. © 2025 Oxford Nanopore Technologies plc. All rights reserved. Oxford Nanopore Technologies products are not intended for use for health assessment or to diagnose, treat, mitigate, cure, or prevent any disease or condition.

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