


Genomic and epidemiological monitoring of yellow fever virus transmission potential

Yellow fever (YF) is an acute viral haemorrhagic disease responsible for up to 60,000 human deaths annually. The disease is caused by a member of the *Flaviviridae* family with transmission primarily through mosquito species that feed on non-human primates. The MinION™ was used in this investigation to characterise genetic diversity and monitor transmission of the virus in real-time, allowing more effective disease surveillance and outbreak control strategies to be implemented¹.

This workflow presented here, as described by Faria *et al*¹ for yellow fever virus (YFV), has been applied to other viruses including Ebola and Zika. More detailed protocol information can be found in Quick *et al.*²

Sample preparation

Sample preparation
~ 9 hours




The Primal Scheme whole genome primer design tool and additional supporting information is available at <http://primal.zibraproject.org>.

Library preparation


Library preparation
~ 1.5 hours

As per Ligation Sequencing 1D (SQK-LSK108) and Native Barcoding (EXP-NBD103) for the MinION with noted amendments*.



*Please note these kits are no longer available. The Ligation Sequencing Kit 1D SQK-LSK108 has since been superseded by kit SQK-LSK109. Two Native Barcoding Expansion kits are now available, providing a total of 24 barcodes (EXP-NBD104 and EXP-NBD114).

Data analysis



In addition to the analysis tools used by Faria *et al*¹, the cloud-based EPI2ME data analysis platform offers a number of real-time analysis pipelines, including species identification and antimicrobial resistance profiling. Find out more at www.nanoporetech.com/analyse.

References

1. Faria, N. R. et al (2018) Genomic and epidemiological monitoring of yellow fever virus transmission potential. *Science*. 361(6405):894-899
2. Quick, J. et al (2017) Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples. *Nat Protoc*. 12(6):1261-1276.
3. GitHub. PoreChop. Available at: <<https://github.com/rwick/Porechop>> [Accessed: 3 May 2018].
4. Li. H. (2013) Aligning sequence reads, clone sequences and assembly contigs with BWA-MEM. arXiv. 1303.3997v1 [q-bio.GN].
5. Loman, N.J. et al (2015). A complete bacterial genome assembled *de novo* using only nanopore sequencing data. *Nature methods*. 12(8):733-5.

Find out more about portable, real-time analysis of viral genomes at www.nanoporetech.com.

*Guppy has superseded Albacore as the production basecaller