

Metagenomic analysis of SARS-CoV-2 respiratory samples via Sequence-Independent Single Primer Amplification (SISPA) and nanopore sequencing.

Samples: SARS-CoV-2 positive nasopharyngeal swabs + negative controls





8

Basecall and demultiplex reads

Perform basecalling and demultiplexing of reads using e.g. MinKNOW™ or guppy



9

Perform host depletion and filtering

Bioinformatically deplete unwanted (e.g. host) reads
Filter reads to discard any under 300 bp in length and/or of Q-score ≤ 7



10

Map reads to SARS-CoV-2 reference and call variants

Map cleaned reads to SARS-CoV-2 reference genome using e.g. minimap2³
Call variants of $\geq 10\%$ frequency using e.g. the callvariants.sh script from the BMap repository⁴
Map reads to reference file containing only coded gene regions using e.g. minimap2³
Call variants in coded gene regions using e.g. the callvariants.sh script from the BMap repository⁴



11

Perform metagenomic analysis

Perform metagenomic analysis of data e.g. via rapid sequence comparison against the GenBank nucleotide database using the sendsketch script from the BMap repository⁴

Kits & devices



Library preparation

Ligation Sequencing Kit (SQK-LSK109)
Native Barcoding Expansion: 1-12 (EXP-NBD104) or 13-24 (EXP-NBD114)



Sequencing

MinION Mk1B or GridION Mk1 + R9.4.1 Flow Cells



References

Sample & library preparation workflow based on protocol described in 2. Bioinformatics workflow based on pipeline described in 5.

- Lewandowski, K. *et al.* (2019). Metagenomic Nanopore sequencing of influenza virus direct from clinical respiratory samples. *J. Clin Microbiol.*, 58:e00963-e01019. doi: 10.1128/JCM.00963-19
- Moreno, G. and O'Connor, D. Protocol: Sequence-Independent, Single-Primer Amplification of RNA viruses V.3. Available at: <https://www.protocols.io/view/sequence-independent-single-primer-amplification-o-bckxiuxn> [accessed 31Mar20]
- Li, H. (2018). Minimap2: pairwise alignment for nucleotide sequences. *Bioinformatics*, 34:3094-3100. doi: 10.1093/bioinformatics/bty191
- BMap repository, GitHub. <https://github.com/BioInfoTools/BMap> [accessed 31Mar20]
- Brown, K. & Moreno, G., SARS-CoV-2 parallel sequencing by Illumina and Oxford Nanopore Technologies, <https://openresearch.labkey.com/Coven/wiki-page.view?name=SARS-CoV-2%20Deep%20Sequencing> [accessed 31Mar20]