

# A complete nanopore-only assembly of an XDR *Mycobacterium tuberculosis* Beijing lineage strain identifies novel genetic variation in repetitive PE/PPE gene regions


With over 10 million new cases reported each year and an increasing incidence of drug resistant infections, tuberculosis (TB) is one of the leading global threats to human health. To better understand the genomic changes that facilitate the emergence and spread of drug resistant *M. tuberculosis* strains (the predominant cause of TB), researchers at the University of Queensland, Australia, utilised nanopore sequencing to deliver a complete *de novo* genome assembly of a highly-transmissible XDR strain from Papua New Guinea<sup>1</sup>.

The MinION™ allowed the full drug-resistance profile to be determined with complete phenotypic concordance. Novel insights were also gained on GC-rich and repetitive genes that are intractable to traditional short-read sequencing technologies.

## Sample preparation

This sample preparation protocol, as described by Bainomugisa *et al*<sup>1</sup> for *M. tuberculosis*, can be applied to other bacterial species.


Sample preparation  
~ 1.5 hours



## Library preparation

As per 1D Genomic DNA sequencing by ligation for the MinION device using SQK-LSK108\* amendments<sup>1</sup>.

Library preparation  
~ 1.5 hours (single sample)  
~ 2 hours (multiplex)



Oxford Nanopore also provides transposase-based library preparation kits, which offer a fast, 10-minute workflow. Both ligation- and transposase-based kits are compatible with barcoding, enabling cost-effective analysis of multiple samples on a single flow cell. For more information, visit [www.nanoporetech.com](http://www.nanoporetech.com).




SEQUENCING


\*SQK-LSK108 has since been superseded by SQK-LSK109.

## Data analysis

### De novo genome assembly



### Variant analysis



In addition to the whole genome assembly and analysis tools used by Bainomugisa *et al*<sup>1</sup>, the cloud-based EPI2ME data analysis platform offers a number of real-time analysis pipelines, including antimicrobial resistance profiling.

### References

- Bainomugisa, A. et al (2018) *Microb Genom.* doi: 10.1099/mgen.0.000188.
- Koren, S. et al (2017) *Genome research.* 27(5):722-36.
- Loman, N.J. et al (2015) *Nature methods.* 12(8):733-5.
- Walker, B.J. et al (2014) *PloS one.* 9(11):e112963.
- Vaser, R. et al (2017) *Genome research.* 27(5):737-746.
- Hunt, M. et al (2015) *Genome biology.* 16:294.
- Kurtz, S. et al (2004) *Genome biology.* 5(2):R12.
- Angiuoli, S.V. et al (2008) *Omics: a journal of integrative biology.* 12(2):137-41.
- Li. H. (2013) arXiv. 1303.3997v1 [q-bio.GN].
- Cingolani, P. et al (2012) *Fly.* 6(2):80-92.
- Kapopoulou, A. et al (2011) *Tuberculosis.* 91(1):8-13.

Find out more about rapid, real-time bacterial sequencing and analysis at [www.nanoporetech.com](http://www.nanoporetech.com).

\*Guppy has superseded Albacore as the production basecaller