

Exploring the potential of adaptive sampling for *Plasmodium* WGS from blood samples


Wim Cuypers^{1*}, Katlijn De Meulenaere^{1,2*}, Julia M. Gauglitz¹, Anna Rosanas-Urgell², Kris Laukens¹, Bart Cuypers¹

1. Adrem Data Lab, Department of Computer Science, Antwerp, Belgium 2. Malaria Unit, Department of Biomedical Sciences, Institute of Tropical Medicine, Antwerp, Belgium | Correspondence: wim.cuypers@uantwerpen.be


Adaptive sampling has the potential to replace more time-consuming protocols for *Plasmodium* WGS

Whole-genome sequencing (WGS) of *Plasmodium* parasites from human blood currently requires **time-consuming lab procedures** to filter out human DNA or enrich *Plasmodium* DNA. We investigated the potential of adaptive sampling to enrich for *Plasmodium* DNA while sequencing unenriched blood samples on a minION device.


Pilot experiment on a dilution series (Kit 12)



Mean 2.7 fold enrichment of *P. falciparum* bases




Faster compared to existing but cheaper methods




Disclaimer

Oxford Nanopore Technologies products are not intended for use for health assessment or to diagnose, treat, mitigate, cure, or prevent any disease or condition.

Direct sequencing of two blood samples (Kit 14)



Competitive read mapping



- ✓ Resistance SNPs -> high concordance with Sanger seq
- ✓ Draft assembly
- ✗ Complete assembly -> 138 & 198 contigs were obtained

Conclusion

Nanopore sequencing in adaptive sampling mode has the potential to replace more time-consuming *Plasmodium*-enrichment protocols and to enable sequencing directly from patient blood, as improvements emerge to increase cost-efficiency.