

Equipment and consumables for Rapid Sequencing Kit family



Oxford Nanopore Technologies
Oxford Science Park, Oxford OX4 4GA, UK
support@nanoporetech.com
www.nanoporetech.com

Equipment and consumables for Rapid Sequencing Kit family

General requirements

Equipment	Microfuge
	Vortex mixer
	Thermal cycler at 30 °C and 80 °C
	P1000 / (P200 / P100) / (P20 / P10) / P2 pipette
	Ice bucket with ice
	Eppendorf tube marker pen
	Timer
Consumables	0.2 ml thin-walled PCR tubes
	Nuclease-free water
	P1000 / (P200 / P100) / (P20 / P10) /P2 pipette tips

Additional requirements by kit

Rapid Low Input by PCR Sequencing Kit (SQK_RLI001)

Third party consumables	Agencourt AMPure XP beads
	LongAmp Taq 2X Master Mix (e.g. NEB M0287)